

CERTIFICATE OF ANALYSIS

For Single Use in Pre-vacuum **STEAM** Sterilizers

This lot of product meets the accepted performance criteria recommended in ISO 11140-1:2014 and ISO 11140-3:2007 for Type 2 Indicators.

LOT 12600

REF A4TS

 : 2020-08







Indicating Ink Origin of Manufacture: United States

 : 2025-08

Indications for Use: 134°C for 3.5 minutes or 121°C for 15 minutes

Requirement	Specification	Result
Color Transition with Effective Vacuum (3.5 Minutes at 134°C or 15 Minutes at 121°C)	Blue to Pink	Pass

Color transition is dependent upon sterilizer size, type, vacuum depth and may vary. A color transition from the initial color blue to a shade of the signal color which is uniformly pink indicates effective steam penetration. If color transition is not uniform or a purple-violet color is achieved, check cycle parameters and equipment function.

Storage and Shelf Life			
	0°C to 40°C		Keep away from sunlight
	20% to 80% Relative Humidity		Keep dry
	Don not re-use; Single use only		Keep away from sterilants; Do not use after expiration date

QUALITY APPROVAL

Date

RETURNS

No returns on dated product will be accepted, except for claim of defect.